

# NEO | SCIENCE

Stay Ahead of the Wellness Curve



## Product Catalogue 2024

(310) 443-8330

[info@neoscience.com](mailto:info@neoscience.com)

[neoscience.com](https://neoscience.com)



# NEO | LIGHT

Go beyond red light with 4 colors and 6 frequencies.  
Create an in-demand and touchless revenue stream.



Boost cellular function



Improve healing



Reduce oxidative stress



Reduce pain



Maintain skin health



Relax muscles



Even skin tone



Improve circulation

## NEO | LIGHT BED Data Sheet

### Size

Weight: 350 kg (771.6 lbs)

Dimension (L × W × H):

92.5 × 41.2 × 39.7 inches

### Product consists of

- **1 Base** (with embedded LEDs, fans and CPU)
- **1 Canopy** (with embedded LEDs, fans, and touch screen controller)
- **1 Variable foot fan** for comfort
- **6ft Power cord**
- **Accessories** including head rest and goggles



### Additional Specifications

#### Electrical Requirements

220-240V, 50/60Hz

30A dedicated circuit breaker

10-gauge wiring

NEMA L6 30R outlet or hard wiring

#### Power consumption

Total rated power : 6336W, 28.8A

Actual power : 3740 W, 21.25A

#### LEDs

31,680 (tripartite spectral diodes)

Pulse option and variable Hz

#### Programs

4 pre-set modes, fully customizable

#### Wavelengths

633 - Red

810 - NIR infrared

850 - NIR infrared

940 - NIR infrared

415 - Blue

532 - Green

#### User Specifications

Max Height : 7' 5"

Max Weight : 661 lb (300 Kg)

# NEO | COOL

100% ELECTRIC LOCALIZED & FACIAL CRYO



The industry is moving away from Nitrogen because it's expensive and difficult to handle. Stay ahead of the curve by switching to our electric localized cryo therapy machine.



## NEO | COOL Data Sheet

### Size

Weight: 55 kg (121.25 lbs)

Dimension (L × W × H):

12.9 × 22.44 × 28.74 inches

### Product consists of

- **1 Base** (with embedded fan and compressor)
- **1 Hose** (with pothook)
- **1 Strut Arm**
- **6ft Power cord**
- **Accessories** (target therapy cup)

### Additional Specifications

#### Electrical Requirements

110 V

7 Amp

#### Power consumption

1000 W

#### Programs

4 pre-set modes, fully customizable

#### Temperature

Lowest Temp -30°C



# BOA Max 2

for business and home

## Lymphatic Drainage Massage Suit for Health, Wellness and Beauty



Immune boosting



Detoxification



Improve circulation



Reduce inflammation



Improve skin health



Relaxation



Supports recovery



Release excess fluids

Great compliment to any facility with body contouring/iv therapy.  
More revenue without more labor.

# BOA Max 2 Data Sheet

### Size

Weight: 13 kg (28.66 lbs)

Dimension (L × W × H):

41 × 34 × 51 inches

### Product consists of

- 1 Suit (with top and bottom)
- 1 Computer System Controller
- 6ft Power cord



### Additional Specifications

#### Electrical Requirements

1110 V

1 Phase Electric, 50/60Hz

1000 W

#### Power consumption

Power supply : 110-230 VAC

#### Programs

Over 20 programs including ipolisis & Cavitation, Sports Massage, Relaxation, Slimming and Cellulite, Lymphatic drainage

#### Pressure

20 up to 140 mmHg

# NEO | LUNA

100% ELECTRIC COLD SAUNA

COMING  
SOON



MADE IN  
USA



## The Luna Difference

- **Seamless Installation**  
Our equipment is designed for effortless setup and relocation.
- **Premium Technology for Nearly Half the Cost**  
Experience cutting-edge performance without the premium price tag.
- **Autonomous Design**  
A completely self-contained unit requiring no external plumbing.
- **Simple 208V Power Requirement**  
Efficiently powered with standard 208V.
- **Proudly Made In The USA and UL Listed**  
Commitment to quality and safety, with 100% American-made craftsmanship.
- **Established Excellence Since 1982**  
A tradition of quality from a manufacturer with over four decades in the industry.
- **24/7 Operational Capability**  
Designed for round-the-clock reliability.
- **Spaciously Designed Interior**  
Comfort meets capacity in a design that's both user-friendly and accommodating.
- **Unmatched Low Operating Costs**  
Optimize your operations with significantly reduced expenses.
- **Continuous Operation with No Interruptions**  
Forget about dry cycles or daily shutdowns.
- **Vortex Cooling Technology**  
Ensures a consistently even distribution of cold for optimal performance.



# NEO | LUNA Data Sheet

## Size

Weight: 680 kg (1500 lbs)

Dimension (L × W × H):

5.5 × 5.5 × 10 feet

## Product consists of

- Insulated Cryogenic chamber
- Refrigeration System
- PLC & touch screen control
- Heated glass door
- Entertainment system

## Additional Specifications

### Electrical Requirements

30-Amp dedicated circuit breaker

125/250 Volt

NEMA L14-30p 4-wire Grounding

Industrial Locking Plug

10-gauge wiring

### Programs

Variable intensity levels

### User Specifications

Max Height :

Max Weight :

Multi-users



# Stand Out With Leading Edge Technology

When you add light therapy, cryo, and lymphatic drainage to your services, you create sustainable, touchless revenue streams for your business.

As wellness business owners ourselves, we help you:

- ✓ Determine which equipment is right for your business or home
- ✓ Get reliable, effective, machine(s) at a competitive price
- ✓ Start using your new product right away with white-glove delivery and easy-to-use preset modes
- ✓ Market your new service to make your new offering profitable as soon as possible



## Position Yourself As An Industry Leader

With any of our health  
optimization products



NEO | LIGHT

NEO | COOL

NEO | LUNA

BOA Max 2

⚙️ Leading Edge Technology

📈 Unlock New Revenue Streams

★ Easy Delivery & Setup

## Choose the Right Equipment

Get the Max Return Quickly



### Jump on a quick call and get a quote

Tell us about your goals, and we'll  
ask a few questions.



### Confidently choose the light bed best suited for you

Get started right away with  
white-glove delivery and simple  
setup.



### Unlock a new revenue stream & grow your business

Provide the red light therapy  
services premier clients are looking  
for.

**Don't wait another week.  
Get the latest technology to attract  
and retain top clients.**

**CALL FOR A QUOTE**



(310) 443-8330

info@neoscience.com

neoscience.com

# NEO | SCIENCE

Stay Ahead of the Wellness Curve